

## Experiment 6

### Thermal expansion of solid bodies

#### 1- Objects of the experiment

- Measuring the linear thermal expansion of a brass tube as a function of the overall length.
- Measuring the linear thermal expansion of glass, steel and brass tubes and determining their linear expansion coefficients.

#### 2- Principles

The length  $s$  of a solid body is linearly dependent on its temperature  $T$ :

$$S = S_0(1 + \alpha.T) \quad (1)$$

$S_0$ : length at  $0^\circ\text{C}$ ,  $T$ : temperature in  $^\circ\text{C}$ .

The linear expansion coefficient  $\alpha$  is determined by the material of the solid body. For a given temperature difference between room temperature  $T_1$  and steam temperature  $T_2$ , the change in length  $\Delta S$  is closely proportional to the overall length  $s_1$  at room temperature:

$$\Delta S \propto S_1 \quad (2)$$

Specifically, we can say:

$$\alpha = \frac{\Delta S}{S_1} \frac{1}{T_2 - T_1} \quad (3)$$

In this experiment, measurements of thermal expansion are conducted on thin tubes through which steam is channeled. The effective length  $S_1$  of each tube can be set as 200, 400 or 600mm by mounting it in a corresponding manner. A dial gauge with 0.01mm scale graduations is used to measure the change in length.

### 3- List of Equipments

| Apparatus                             | Catalogue Number |
|---------------------------------------|------------------|
| 1 Expansion apparatus                 | 38134            |
| 1 Holder for dial gauge               | 38136            |
| 1 Dial gauge                          | 36115            |
| 1 Steam generator, 550 W/230 V        | 30328            |
| 1 Silicone tubing, dia. 7X1.5 mm, 1 m | 667194           |
| 1 Petri dish, 150 X 25 mm             | 664185           |
| 1 Thermometer -10°C to 110°C          | 38234            |

### 4- Setup

- Set up the experiment as shown in Fig. 1.




Fig. 1 Experiment setup for measuring the linear thermal expansion with the expansion apparatus

- Screw on the holder for dial gauge (e) and clamp the dial gauge in place.
- Attach the fixed bearing (a) of the expansion apparatus at the 600 mark and slide the open end of the brass tube into the fixed bearing.

- Slide the closed end of the brass tube into the guide fitting (b) so that the hose nipple (f) is pointing laterally downwards.
- Tighten the screw to fix the brass tube in the fixed bearing (the screw must engage the ring groove of the tube).
- Insert the extension piece (c).
- Cut off a 20cm length of silicone tubing, slide the section over the hose nipple (f) and place a Petri dish below this to catch the condensation.
- Use the long tubing section to connect the open end of the brass tube to the steam generator.

### **5- Carrying out the experiment**

- Determine the room temperature  $T_1$  and write this down.
- Read off and write down the zero position of the dial gauge.
- Fill the steam generator with about 2cm of water, close the apparatus and plug it in.
- Read off and write down the maximum pointer deflection of the dial gauge.
- Allow the brass tube to cool down to room temperature.
- Attach the fixed bearing of the expansion apparatus at the 400 mark and tighten the screw so that it engages the ring groove of the tube.
- Refill the steam generator with water, check the zero position of the dial gauge and repeat the measurement.
- Move the fixed bearing to the 200 mark and repeat the experiment.
- Replace the brass tube with the steel tube; attach the fixed bearing at the 600 mark and repeat the measurement.
- Conduct the same measurement on the glass tube.
- Complete the following table.

Table 1: Linear expansion  $\Delta S$  between room temperature  $T_1$  and steam temperature  $T_2$  as a function of the effective length  $S_1$  of the material.

| Material | $S_1$ (mm) | $\Delta S$ (mm) |
|----------|------------|-----------------|
| Brass    | 600        | -----           |
| Brass    | 400        | -----           |
| Brass    | 200        | -----           |
| Glass    | 600        | -----           |
| Steel    | 600        | -----           |

- Determine  $T_1$  ,  $T_2$  and  $T_2 - T_1$ .
- Prepare a sheet of graph paper for plotting  $\Delta S$  of a brass tube versus  $S_1$ . You should make  $\Delta S$  the vertical axis and  $S_1$  the horizontal axis. Each axis should be labeled and appropriate units indicated. The graph should have a title.
- Plot your data on the graph.
- Determine the slope and the y-intercept point by using the least square method (see appendix).
- Draw best fit line to the points on your graph.
- Use the measuring results of table 1 and equation (3) to determine the linear expansion  $\alpha$  for different materials.

Table 2. Linear expansion coefficient  $\alpha$  for different materials

|          |            | Measurement           | Literature            |
|----------|------------|-----------------------|-----------------------|
| Material | $S_1$ (mm) | $\alpha$ ( $K^{-1}$ ) | $\alpha$ ( $K^{-1}$ ) |
| Brass    | 600        | -----                 | $18 \cdot 10^{-6}$    |
| Steel    | 600        | -----                 | $11 \cdot 10^{-6}$    |

|       |     |       |                   |
|-------|-----|-------|-------------------|
| Glass | 600 | ----- | $3 \cdot 10^{-6}$ |
|-------|-----|-------|-------------------|

- Compare these results with literature values and discuss your results.